



” 10 cm (4 inches) - that’s the difference between falling or rolling out of bed. Such a small distance, yet these 10cm can avoid extended hospital stays and direct costs.

## FloorLine-i

Floor-level-bed for use in hospitals and long-term care facilities

- Minimum height from the floor of only 10 cm
- Can be used both in hospitals and long term care facilities
- Minimizes patient restraints and use of side rails

The FloorLine-i is a sound purchase for both hospitals and long term care facilities. It features standard bed positions such as Trendelenburg and cardiac chair position providing excellent care and can be equipped with a variety of optional accessories to accommodate multiple patient requirements.

The FloorLine-i is an optimal floor-level-bed solution that promotes both patient as well as caregiver safety. Its minimum height from the floor is just 10 cm (3.9”) and at the same time no compromises are made for caregivers with its ergonomic 80cm (31.5”) working height.



Part No. (including hand control and end boards)	Platform Position min- max	Sleeping Dimensions L x W	Bed Dimensions L x W	Functions	Compliance Standards	Maximum User Weight/ Safe Working Load	Product Weight
Europe: 90300 Australia : 90301 North America: 90302 UK: 90303	10cm - 80cm (4" - 31.5")	2000mm x 900mm (78.7" x 35.5")  With Extension Kit 2175mm x 900mm (85.5" x 35.5")	2325mm x 925mm (91.5" x 36.5")  With Extension Kit 2500mm x 925mm (98.5" x 36.5")	Twin castors Battery Back Up, Auto Contour Button for Cardiac Chair Position, Hand- set w Lockout	IEC60601-2-52:2009 EN1970:2002 AS/NZS 3200.2.38:2007	Maximum user weight: 185kg (407lbs)  Safe working load: 250kg (550lbs)	138 kg (303 lbs)

FloorLine-i can be fitted with a variety of optional accessories such as dedicated side-rails, self-help pole and length extension kits.