



“ The FloorLine LTC with its home environment design is the optimal choice for long-term care facilities

FloorLine LTC

Floor-level-bed for long term care facilities

- Lowest height from the floor of only 9.5 cm (3.7”)
- The optimal choice for long term care
- Easy to fold, store and transport

The FloorLine LTC with its home environment design is the optimal choice for long term care facilities. Its one-piece design makes it easy to fold and store away efficiently. The FloorLine LTC features all the standard required positions often found in hospital beds such as Trendelenburg, reverse Trendelenburg and cardiac positioning.

The FloorLine LTC is the perfect solution for preventing fall injuries. Its minimum height from the floor is just 9.5 cm(3.7”) providing a safe and considerate way of taking care of people who may be at risk of falling out of their bed.



Part No.	Platform Position min- max	Sleeping Dimensions L x W	Bed Dimensions L x W	Functions	Compliance Standards	Maximum User Weight/ Safe Working Load	Product Weight
Europe: 90600 UK: 90602 North America: 90604 USA class 1: 90606 Australia: 90608 Europe with wooden end boards: 90601 UK with wooden end boards: 90603 NA with wooden end boards: 90605 NA class 1 with wooden end boards: 90607 AUS with wooden end boards: 90609	9.5 cm - 70 cm (3.7” - 27.5”)	2000mm x 900mm (78.7” x 35.5”) With Extension Kit 2175mm x 900mm (85.5” x 35.5”)	2150mm x 905mm (84.5” x 35.5”) With Extension Kit 2325mm x 905mm (91.5” x 35.5”)	Auto Contour Button for Cardiac Chair Position, Handset w Lockout (Batteribackup as accessory, part no: 75805)	IEC 60601-2-52:2009 EN1970: 2002 AS/NZS 3200.2.38: 2007	Maximum user weight: 150kg (330lbs) Safe working load: 200kg (440lbs)	106 kg (234 lbs)

FloorLine LTC can be fitted with a variety of optional accessories such as dedicated side-rails, self-help pole, back up battery and length extension kits.