

Padded Transfer Commode Chair



Style# MHPTCM



Padded Transfer Commode Chair: MHPTCM

- o Dual-use Transfer bench and Commode Seat
- o New design with easy access for cleaning
- o Comfortable Cushioned Seat and Backrest
- o Aluminum Construction that is Durable, Lightweight, and Corrosion Free
- o Suction Cup Feet for added stability
- o Commode Bucket included
- o Weight capacity 300lbs

MHPTCM



MOBB Health Care Ltd.

116A Viceroy Rd. Concord,
ON L4K 2M1 Canada
Tel: 905-669-0348
Fax: 905-669-2796
www.mobbhhc.com

SPECIFICATION

Weight Capacity	300 lbs
Warranty	Limited One Year
Material	Aluminum

Measurements

