

## Gravity 6 Long Term Care Pressure Redistribution Mattress

HCPCS: E0184,E0184



- High quality latex free therapeutic foam mattress provides pressure redistribution and shear/friction reduction
- Deluxe horizontal, cross-cut foam mattress provides comfort, support, and pressure redistribution over 5 therapeutic pressure zones
- Better compliance and durability than innerspring mattresses allows for improved patient satisfaction and fewer service calls
- Non-skid waterproof vinyl bottom cover keeps mattress safely in place
- Heel section is sloped offloading pressure from vulnerable heels
- Removable multi-stretch cover is fluid resistant, low shear and vapor permeable helping protect the skin from friction and moisture
- Durable cover tolerates disinfectant cleaning supporting infection control, refer to user manual for suggested disinfectants
- Meets flammability standard CFR 1632 and CFR 1633

PRODUCTS	Description
15670	GRAVITY 6 MATTRESS, 80", 1/EA
SPECIFICATIONS	Description
Warranty	3 Year Foam, 3 Year Cover
Weight Capacity	350 lbs