

## Prime Plus Care Bed Model P750

HCPCS:



- 750 lbs weight capacity
- 700 lbs Net weight capacity
- 11" low height 30" high working height
- Push button integrated width expansion from 39" to 48" width
- Integrated 80", 84" and 88" sleep surface length helps comfortably accommodate the tallest residents
- One touch floor lock
- Quiet, smooth operation/Easy-to-use Hand Pendant
- Cardiac Comfort Chair
- 4 heavy duty locking casters roll easy at any height and lock secure
- Fixed adjustable wall bumper
- Warranty: 15 years on frame lifetime on welds 5 years on electronics
- Options:
- Deluxe Assist Rails, Assist Bars
- Bed ends available in 2 styles and 3 colors
- Trapeze
- Battery Backup
- Bed Transporter

PRODUCTS	Description
P750	LTC Bariatric Bed 750 lbs

  

SPECIFICATIONS	Description
Warranty	15 Year on Frame, 5 Year on Electrical and Mechanical. Lifetime on All Welds
Weight Capacity	750 lbs
Built-in Width Expansion	39", 42", 48"
One Step Locking	Yes
Sleep Surface	80" - 84" - 88"

  

ACCESSORIES	Description
P903-T-BBK	BATTERY BACKUP W/MOUNTING, P903, Timotio, EA
PLTCAB	LTC Assist Bar, 1PR, PR
PLTCAR	LTC Assist Rail, 1PR, PR

  

PARTS	Description
SP01-P903-10006E	Head Board Bracket Knobs , P903, 1/pr, PR
SP01-P750-10004E	P750 head deck, 1/EA, EA
SP01-P750-10005E	P750 thigh deck, 1/EA, EA

**PARTS****Description**

SP01-P750-10006E	P750 Foot Deck, 1/EA, EA
SP01-PO10040209E-G	Brake Pedal, Green, P903, 1/ea, EA
SP01-PO10040209E-R	Brake Pedal, Red, P903, 1/ea, EA
SP01-P903-10008E	Head/FootBoardMountingScrews,P903,1/pkg, PKG
SP01-P703-60-02E	Power Cord,P503/703 & P903, EA
SP01-TH12-1060-001	Hand Pendant, P903, 1/ea W/ backlit, EA