


# FUSION<sup>2K</sup>™

Featuring SelectAir® Therapy

**Wicking Action:** Air escapes through tiny holes on the top of the mattress and circulates under a vapor-permeable cover to wick moisture away from patient.



## iSense Processing™

Advanced microprocessors continually monitor patient's immersion.

## FUSION<sup>2K</sup>™ Features

### 'True' Low Air Loss Technology

Allows air to escape from the top of the mattress and directly interact with moisture at the patient/bed boundary.

### Pulsation

Provides two cycle settings that continually readjust the system to ensure patient comfort.

### Firm Perimeter Support

For ease of transfer and safety.

### User-Friendly Control Unit

Rapid inflation, easy set-up and a motor so quiet, the patient will hardly know it's there!

**Zoned, "true" low air loss used in acute care hospitals nationwide with firm perimeter support for ease of transfer and safety**

**Weight Limit:** 300 lbs. (136 kg)

**Mattress Size:** 80" (or 84") x 36" x 6"

**Cover:** Water resistant, vapor-permeable, flame resistant, biocompatible and quilted

**Warranty:** Two Year Limited Warranty on control unit and soft goods

LIT-0007 | R100121



**Blake**  
Medical Distribution

42 Niagara Street,  
Hamilton, ON L8L 6A2  
[www.blakemedical.ca](http://www.blakemedical.ca)

Phone: 1 (866) 662-3459  
Fax: 1 (905) 560-1342  
Email: [blake.orders@blakemedical.ca](mailto:blake.orders@blakemedical.ca)