

# SELECTAIR® MAX

Low Air Loss Mattress Replacement



## User-Friendly Control Unit

Rapid inflation, easy set-up and a motor so quiet, the patient will hardly know it's there!

## The SelectAir® System

Quietly circulates air through individual cells and allows air to escape from the mattress.

**Wicking Action:** Air escapes through tiny holes on the top of the mattress and circulates under a vapor-permeable cover to wick moisture away from patient.

## SELECTAIR® MAX Features

**Pulsation (alternating pressure):** Provides three cycle settings that continually readjust the system to ensure patient comfort.

**System Lock:** Helps to prevent accidental adjustments and automatically resets the system in case of a power failure.

**Auto Comfort:** Takes the guess work out of programming. During set-up, mattress sets optimal pressure according to patient's weight distribution.

**Zoned, "true" low air loss used in acute care hospitals nationwide**

**Weight Limit:** 500 lbs. (227 kg)

**Mattress Size:** Fits most standard bed frames. Custom sizes available. Inquire for details.

**Cover:** Water resistant, vapor-permeable, flame resistant and biocompatible

**Warranty:** Two Year Limited Warranty on control unit and soft goods

LIT-0003 | R100121



**Blake**  
Medical Distribution

42 Niagara Street,  
Hamilton, ON L8L 6A2  
[www.blakemedical.ca](http://www.blakemedical.ca)

Phone: 1 (866) 662-3459  
Fax: 1 (905) 560-1342  
Email: [blake.orders@blakemedical.ca](mailto:blake.orders@blakemedical.ca)