


Adapt 2 Paediatric

Modular system



At a glance

- Suitable for young children with medium to high risk for pressure injuries
- Modifiable surface allows for addition of up to three ROHO neonatal isolette inserts
- Fits exactly in a standard crib
- Full enclosure cover with Recovery 5 top
- Maximum client weight of up to 200 lb
- 30"W x 80"L x 6"H
- 24-month warranty

The Adapt 2 Paediatric mattress with optional ROHO® isolette or Serene foam inserts is ideal for vulnerable children with high to very high risk for pressure injuries. The surface inserts can be configured based on the individual client's needs (see samples ideas below) and can hold up to three (3) ROHO isolette inserts. This modified mattress fits exactly in a standard crib. The 1.5" wide firm perimeter borders and 4" permanent gel foam foot border hold the inserts in place without the need for snaps. The mattress base is protected by a Recovery 5 anti-microbial incontinent cover, as well as a Recovery 5 protective top cover.

