

# Align and Align Bariatric Mattress

The Align surface with SRT comes with firm borders, a raised sleep surface, and full enclosure cover with a Recovery 5 top. It is a unique multi-zone surface that provides differing immersion and envelopment levels that match the contours of the body. There is one piece of HR foam at the torso, with 4" square cuts in the foam under the sacrum, and 2" square cuts under the heels to release the surface tension. The foam density and weight capacity are what set our Align & Align Bariatric options apart. Upgrade the mattress with a raised bolster or R5 Ultra cover.

**Note:** Custom sizes available, please call for details.

## Features at a glance

- Ideal for long-term care clients
- Multi-zone surface with differing immersion and envelopment levels
- Wide transfer edges to help reduce the risk of falls
- Shear Reduction Technology (SRT)
- Full enclosure cover with Recovery 5 top
- Maximum client weight of up to 350-500 lb for regular Align, and 600-800 lb for bariatric
- Available sizes:
  - 36" / 39" / 42" / 48" / 54" / 60" x 80"L x 7"H
  - 36" / 39" / 42" x 84"L x 7"H
  - Custom sizes also available
- 24-month warranty

## This mattress is

