

# Atlas Auto Comfort Mattress

Atlas Auto Comfort uses equalization technology to provide maximum pressure redistribution to prevent pressure injuries or to support the healing process. Beneath the 1" top layer of convoluted high resiliency foam there are interconnected chambers in a closed circuit to provide self-adjusting pressure redistribution. Clients also benefit from 4" wide firm borders to ensure a solid transfer edge. Atlas Auto Comfort includes a 5-degree sloped foot section to help transfer excess pressure from the heels.

## Features at a glance

- Suitable for vulnerable clients at risk for pressure injuries
- Interconnected chambers in a closed circuit allow for self-adjusting pressure redistribution
- 5-Degree Sloped foot section helps reduce pressure from heels
- 1" high convoluted foam top layer for comfort
- 4" wide firm borders for safe transfers
- Fluid-resistant, four-way stretch cover
- Supports clients weighing up to 500-700 lb., dependent on size of surface
- Available sizes:
  - 35" / 39" / 42" x 80"L x 7"H
  - Custom sizes also available

