

Atlas Hybrid

The Atlas Hybrid is a non-powered, self-adjusting mattress that utilizes equalization technology to provide maximum pressure redistribution for the prevention and treatment of pressure injuries. The multi-zone top has Serene™ foam in the middle section for immersion plus

a 5-degree sloped foot section made with Super Soft™ foam to help transfer excess pressure from the heels. Beneath the top layer, interconnected chambers provide pressure redistribution, and there are 4" wide firm borders to ensure a solid transfer edge.

A manual ROHO® air pump can be used to easily adjust the softness/firmness of the chambers. Both the closure kit and pump come included with the Atlas Hybrid.

Features at a glance

- Suitable for vulnerable clients at risk for pressure injuries
- Interconnected chambers and valves for pressure redistribution
- ROHO pump adjusts the firmness
- Sloped foot section helps reduce pressure from heels
- 4" wide firm borders for safe transfers
- Fluid-resistant, four-way stretch cover that is wipeable for easy cleaning
- Maximum client weight of up to 500-700 lb depending on size of surface
- Available sizes:
 - 35" / 39" / 42" x 80"L

