

Therapeutic Surface

Domus Auto

Alternating Pressure Redistribution System

The Domus Auto operates on a 10 minute cycle which allows the bladders to inflate/deflate allowing for optimal offloading and transfer of body pressure. The system features Alternating Pressure Therapy and Constant Low Pressure (static), along with key features such as max inflate and seat inflation for the Fowler position and an automatic adjustment system allowing for fast and simple set-up. For added protection the Domus Auto contains multiple safety features such as cell-on-cell bladder to avoid bottoming out, a highly visible CPR release valve, and an auto-lock on the settings after five (5) minutes

Features at a glance

- Optimal for clients with medium to high risk for pressure injuries
- Automatic adjustment system for fast and simple setup
- Low Air Loss for microclimate control
- Lower 3" of air cells for the whole surface remain static to prevent bottoming out
- Silver+ nylon, fluid-resistant, quilted cover
- Zippered top cover allows for one-step removal for easy cleaning
- Visual and audible alarm enhances clients' safety & security
- Supports clients with a weight of up to 450 lb
- **Optional features:** Add-on bolsters with full enclosure R5 covers and 6" wide firm transfer edges transform Domus Auto systems into a 48" wide Domus Edge surface to meet your clients' needs
- Available sizes:
 - 35.5" W x 79" L x 8"H
- 36-month warranty

