

Therapeutic Surface

Domus 4

Domus 4 Alternating Pressure Therapy operates on a 10 minute cycle which allows the bladders to inflate/deflate allowing for optimal offloading and transfer of body pressure. The system features Alternating Pressure Therapy and Constant Low Pressure (static), along with key features such as max inflate and seat inflation for the Fowler position. For added protection the Domus 4 contains multiple safety features such as cell-on-cell bladder to avoid bottoming out, a highly visible CPR release valve, and an auto-lock on the settings after five (5) minutes.

Features at a glance

- Optimal for clients with medium to high risk for pressure injuries
- Low Air Loss for microclimate control
- Cell-on-cell bladder to avoid bottoming out
- Silver+ nylon, fluid-resistant, quilted cover
- Zippered cover allows for one-step removal
- Visual and audible alarm enhances clients' safety & security
- Maximum client weight of up to 450 lb
- **Optional features:** Add-on bolsters with full enclosure R5 covers and 6" wide firm transfer edges transform Domus 4 systems into a 48" wide Domus Edge surface to meet your clients' needs.
- Available sizes:
 - 35.5" W x 79" L x 8"H
 - Gap filler **extends surface to 84" L**
- 36-month warranty

