

## Walker Tray with Cup Holders



HCPCS:



FIGURE A



FIGURE B

- Allows personal items to be carried from room to room
- Contains two cup holders
- Fits most manufacturers' 2-button walkers
- Easily folds down when not in use
- Made of easy-to-clean, durable plastic
- Easy-to-install
- Size: 16" (L) x 12" (W) x 2" (H)

**PRODUCTS**

**Description**

10125 Walker Tray, 1/bx, RTL

**SPECIFICATIONS**

**Description**

Warranty Limited 1 Year

Dimensions 16" (L) x 12" (W) x 2" (H)