

Winnie Lite Supreme/Go Lite 3 Wheel Aluminum Rollator

HCPCS: E0143,E0143




FIGURE A




FIGURE B



FIGURE C



FIGURE D

- Special loop lock made of internal aluminum casting operates easily and ensures safety (Figure A)
- 7.5" casters with soft-grip tires are ideal for indoor and outdoor use
- Comes standard with carry pouch
- Handles easily adjust in height with self-threading knob (Figure B)
- Easy, one-hand folding

PRODUCTS

Description

RTL199	Winnie Lite Supr.w/TanPlaidPou,1/cs
--------	-------------------------------------

SPECIFICATIONS

Description

Warranty	Limited Lifetime
Handle Height	30.75" - 37.5"
Weight Capacity	300 lb
Carton Shipping Weight	15 lb
Weight	11 lb
Wheel Size	7.5"

PARTS

Description

10208KN	Knob,Ht.Adj.Rollator,w/Attch.1/ea, EA
9500001232	Brake Cable/Adjust,Lower,1/ea, EA
9500004101	Cable,Fastner,1/ea, EA
950330049901	Hand Grip, Right, 3-wheel Rollator, 1/ea, EA
950330039901	Hand Grip, Left, 3-wheel Rollator, 1/ea, EA
950000199901	Generic Brake Cable Only for Rollat,1/ea, EA