

Unity

# Fall Mat

The Unity Fall Mat is a strong and resilient injury mitigation mat designed to provide protection beside the bed and reduce the impact during a fall event.



## Features at a glance

- **Full bed length protection:** 72" length by 24" wide
- **Strong and Resilient:** Polyurethane foam inner core to provides maximum shock absorption.
- **Low profile with tapered edges:** Does not impede access to the bed and minimizes the tripping hazard.
- **Anti-slip bottom:** Does not move on tiles, carpet, or vinyl flooring.
- **Textured top:** Increases surface traction
- **Easy to clean:** Can be wiped with soap and water, or hospital grade cleaners.