

Permobil

# Heavy duty products



permobil





Up to 8" of  
seat elevation

Base width: 26"  
Base length: 38"

# Permobil M300 Corpus® HD

The M300 Corpus HD takes all the things you love about Permobil’s mid-wheel drive and amplifies them. Outfitted with stronger actuators, more powerful motors, and a reinforced Corpus seating system, it’s built to withstand the rigors of an active heavy duty lifestyle. Designed specifically for users up to 450 lb, the suspension has been strengthened to provide you with more stability and completed with elegant features like fully adjustable legrests and up to four power seat functions.



Weight capacity— 450 lb

## Specifications

Max user weight	450 lb
Max speed	5 mph
Range	13 mi (Grp 34)   20 mi (Grp 24)
Base width	26"
Base length	38"
Seat-to-floor height	18", 19", 20" (elevator) 18", 19", 20" (fixed)
Base turning radius	22"
Ground clearance	3"
Recommended battery type	Group 34 or 24
Weight including batteries	425 lb (Group 34) 447 lb (Group 24)
Drive electronics	R-Net 120A

Optional seat elevator	8"
Tilt options	Posterior: 0° – 45°
Backrest angle - power recline	85° – 150°
Legrest angle - power	85° – 175°
Seat depth	19" – 25" (by 1" incr.)
Seat width	19", 21", 23", 25"
Backrest height	21", 24" – 29" (by 1" incr.)

*\* Actual driving range will vary based on driving and battery conditions.  
All specifications may vary based upon configuration.  
Chair weights are approximate and include base, batteries, power seat functions and seating system without accessories.*

## Color options

Volcano red	Pop Star pink	Sublime green	Cobalt blue	Purple rain	Platinum silver	Phantom grey	Onyx black



# Permobil M3 Corpus

Meticulously engineered to help you maintain better positioning, feel more secure and reduce driving fatigue so you can conquer daily activities with ease. Its compact base offers maneuverability in tight spaces without sacrificing performance. Featuring the innovative FlexLink suspension geometry with independent caster swing-arms helping to keep you stable and in control as you traverse the obstacles in your world at up to 6 mph.




## Specifications

Max user weight	330 lb
Max speed	6 mph
Range	16 mi (Grp 34)   20 mi (Grp 24)
Base width	24"
Base length	35.5"
Seat-to-floor height	17.5", 18.5", 19.5"
Base turning radius	20.5"
Ground clearance	3"
Recommended battery type	Group 34 or 24
Weight including batteries	385 lb (Group 34) 405 lb (Group 24)
Drive electronics	R-Net 120A
Optional seat elevator	12"


Tilt options	Posterior: 0° – 50° ActiveReach: 5°, 10° & 20° (up to 300 lb)
Backrest angle - power recline	85° – 180°
Legrest angle - power	85° – 180°
Seat depth	14" – 22" (1" incr.)
Seat width	17" – 23" (2" incr.)
Backrest height	20", 23" – 28" (1" incr.)

*\* Actual driving range will vary based on driving and battery conditions.  
All specifications may vary based upon configuration.  
Chair weights are approximate and include base, batteries, power seat functions and seating system without accessories.*

## Color options



## Premium options





# Permobil M5 Corpus

Take control with high-performance motors enabling a top speed of 7.5 mph and exceptional low-end torque that lets you pivot like a pro. Experience unsurpassed stability and ride comfort with FlexLink Pro independent drive wheel suspension. The M5 Corpus will take you to new heights with increased stability and full suspension in seat elevation.

Up to 14" of ActiveHeight™ & 20° of ActiveReach™

Base width: 25"  
Base length: 37.5"



Weight capacity— 330 lb

## Specifications

Max user weight	330 lb
Max speed	7.5 mph
Range	20 mi
Base width	25"
Base length	37.5"
Seat-to-floor height	17.5", 18.5", 19.5"
Base turning radius	21.5"
Ground clearance	3"
Recommended battery type	Group 24
Weight including batteries	423 lb
Drive electronics	R-Net 120A
Optional seat elevator	14"

Tilt options	Posterior: 0° – 50° ActiveReach: 5°, 10° & 20° (up to 300 lb)
Backrest angle - power recline	85° – 180°
Legrest angle - power	85° – 180°
Seat depth	14" – 22" (1" incr.)
Seat width	17" – 23" (2" incr.)
Backrest height	20", 23" – 28" (1" incr.)

\* Actual driving range will vary based on driving and battery conditions.  
All specifications may vary based upon configuration.  
Chair weights are approximate and include base, batteries, power seat functions and seating system without accessories.

## Color options

Chili

Sunset

Canary

Spring

Aqua

Sky

Violet

Flamingo

Arctic

Midnight

Gold

Platinum

Bronze

## Premium options



# Permobil F3 Corpus

The F3 Corpus with a 6 mph top-end speed features optimized battery weight distribution, lower center of mass, and fully adjustable ComfortRide suspension that improve stability while traversing inclines or declines. Experience the enhanced functionality that combines outstanding chassis performance with Permobil's industry-leading Corpus seating system.



## Specifications

Max user weight	330 lb	Drive electronics	R-Net 120A
Max speed	6 mph 5 mph <i>without support wheels</i>	Optional seat elevator	12"
Range	16 mi (Grp 34)   20 mi (Grp 24)	Tilt options	Posterior: 0° – 50° ActiveReach: 5°, 10°, 20° & 30°
Base width	24"	Backrest angle - power recline	85° – 180°
Base length	40" with support wheels 36" without support wheels	Legrest angle - power	85° – 180°
Seat-to-floor height	17.5", 18.5", 19.5"	Seat depth	14" – 22" (1" incr.)
Base turning radius	26.5"	Seat width	17" – 23" (2" incr.)
Ground clearance	3"	Backrest height	20", 23" – 28" (1" incr.)
Recommended battery type	Group 34 or 24	<i>* Actual driving range will vary based on driving and battery conditions. All specifications may vary based upon configuration. Chair weights are approximate and include base, batteries, power seat functions and seating system without accessories.</i>	
Weight including batteries	386 lb (Group 34) 407 lb (Group 24)		

## Color options

Chili

Sunset

Canary

Spring

Aqua

Sky

Violet

Flamingo

Arctic

Midnight

Premium options

Gold


Platinum

Bronze



# Permobil F5 Corpus

Power meets performance on the F5 Corpus with a 7.5 mph top-end speed and upgraded ComfortRide Pro suspension. All-wheel independent, fully adjustable, oil dampened shocks provide the ultimate in ride comfort and reduced vibration while driving in all seat positions, including an unprecedented 14” ActiveHeight seat elevation. Climb and reach higher with increased obstacle negotiation, exceptional motor torque, and enhanced stability allowing up to 45° of ActiveReach forward tilt functionality.



## Specifications

Max user weight	330 lb
Max speed	7.5 mph 6 mph without support wheels
Range	20 mi
Base width	25.5"
Base length	42" with support wheels 38" without support wheels
Seat-to-floor height	17.5", 18.5", 19.5"
Base turning radius	30"
Ground clearance	3"
Recommended battery type	Group 24
Weight including batteries	412 lb
Drive electronics	R-Net 120A

Optional seat elevator	14"
Tilt options	Posterior: 0° – 50° ActiveReach: 5°, 10°, 20°, & 45°
Backrest angle - power recline	85° – 180°
Legrest angle - power	85° – 180°
Seat depth	14" – 22" (1" incr.)
Seat width	17" – 23" (2" incr.)
Backrest height	20", 23" – 28" (1" incr.)

*\* Actual driving range will vary based on driving and battery conditions.  
All specifications may vary based upon configuration.  
Chair weights are approximate and include base, batteries, power seat functions and seating system without accessories.*

## Color options



## Premium options



# TiLite Aero T

Our most advanced aluminum chair is uniquely modeled and handmade to each user's specifications. The dual-tube frame of the Aero T provides superior rigidity, while minimizing weight and maximizing performance. Thanks to the TiFit process, you get an optimized fit plus adjustability that meets your needs.



Note: The ROHO® AGILITY® Active Carbon back supports up to 265 lb.



















## Specifications

Max user weight	265 lb standard 350 lb heavy duty	Camber	0°, 2°, 4°, or 6°
Frame features	Titanium dual-tube rigid, TiFit, full adjustability	Center of gravity adjustment	5"
Total chair width*	23"	Front seat width (seat taper)	0", 1", 2", 3"
Seat width	12" – 20"	Precision lock folding seat back	Adj. height: 9.5" – 22" Adj. angle: 80° – 101°
Seat depth	12" – 20"	Ergonomic seat	6", 7", 8"
Front seat height	15" – 21"	Transit options	Q-Straint lockdown approved
Rear seat height	12.5" – 21"	Transport weight*	12 lb
Front frame angle	70°, 75°, 80°, 85°, or 90°	Warranty	Five years on frame One year on TiLite components








\*Weights and widths based on 16" x 16" chair with 0° camber, minimum rear wheel spacing and all standard features. Rear wheels are not included in transport weight. For reference, a pair of standard 24" rear wheels weighs 7.2 lb.

## Frame colors

High gloss paint			Candy paint		Metallic paint			
								
TiLite red	Neon pink	High gloss black	Ruby red	Azure blue	Electric plum	Tangerine metallic	Acid green metallic	Ocean blue metallic
Matte textured paint			Matte paint		Pearlescent paint			
								
Flame orange	OD green	Hyper blue	Titanium grey	Black lava	Midnight blue pearl	White pearl		

## Component colors

# Anodized

		
Red	Orange	Gold
		
Green	Blue	Purple
		
Black		





# TiLite TRA

With the power and style of its titanium dual-tube frame, plus TiFit precision for an exact fit, the TRA provides performance and adjustability made just for you.





















## Specifications

Max user weight	265 lb standard 350 lb heavy duty	Camber	0°, 2°, 4°, or 6°
Frame features	Titanium dual-tube rigid, TiFit, full adjustability	Center of gravity adjustment	5"
Total chair width*	23"	Front seat width (seat taper)	0", 1", 2", 3"
Seat width	12" – 20"	Precision lock folding seat back	Adj. height: 9.5" – 22" Adj. angle: 80° – 101°
Seat depth	12" – 20"	Ergonomic seat	6", 7", 8"
Front seat height	15" – 21"	Transit options	Q-Strait lockdown approved
Rear seat height	12.5" – 21"	Transport weight*	11.1 lb
Front frame angle	70°, 75°, 80°, 85°, or 90°	Warranty	Five years on frame One year on TiLite components








\*Weights and widths based on 16" x 16" chair with 0° camber, minimum rear wheel spacing and all standard features. Rear wheels are not included in transport weight. For reference, a pair of standard 24" rear wheels weighs 7.2 lb.

## Frame colors

High gloss paint			Candy paint		Metallic paint			
								
TiLite red	Neon pink	High gloss black	Ruby red	Azure blue	Electric plum	Tangerine metallic	Acid green metallic	Ocean blue metallic
Matte textured paint			Matte paint		Pearlescent paint		Natural titanium	
								
Flame orange	OD green	Hyper blue	Titanium grey	Black lava	Midnight blue pearl	White pearl	Satin	Polished

## Component colors

# Anodized

		
Red	Orange	Gold
		
Green	Blue	Purple
		
Black		





# TiLite ZRA

With full adjustability and a minimalist mono-tube frame, the ZRA dishes out style without boundaries. With its TiFit titanium frame, you have the complete package: strength, style, and the ability to set up and configure your chair just for you.

Weight capacity— 300 lb





















## Specifications

Max user weight	265 lb standard 300 lb heavy duty	Camber	0°, 2°, 4°, or 6°
Frame features	Titanium mono-tube rigid, TiFit, full adjustability	Center of gravity adjustment	6"
Total chair width*	23"	Front seat width (seat taper)	0", 1.5", 2.5"
Seat width	10" – 20"	Precision Lock folding seat back	Adj. height: 9.5" – 22" Adj. angle: 80° – 101°
Seat depth	10" – 20"	Ergonomic seat	6", 7", 8"
Front seat height	16" – 21"	Transit options	Transit tie-down
Rear seat height	12.5" – 20"	Transport weight*	10.7 lb
Front frame angle	70°, 75°, 80°, 85°, or 90°	Warranty	Five years on frame One year on TiLite components








\*Weights and widths based on 16" x 16" chair with 0° camber, minimum rear wheel spacing and all standard features. Rear wheels are not included in transport weight. For reference, a pair of standard 24" rear wheels weighs 7.2 lb.

## Frame colors

High gloss paint	Candy paint	Metallic paint
 TiLite red	 Ruby red	 Tangerine metallic
 Neon pink	 Azure blue	 Acid green metallic
 High gloss black	 Electric plum	 Ocean blue metallic
Matte textured paint	Matte paint	Pearlescent paint
 Flame orange	 Titanium grey	 White pearl
 OD green	 Black lava	 Satin
 Hyper blue	 Midnight blue pearl	 Polished

## Component colors

# Anodized

		
Red	Orange	Gold
		
Green	Blue	Purple
		
Black		





# TiLite ZR

With the ZR’s minimalist mono-tube frame, less truly is more. The clean aesthetic of the ZR turns heads while the titanium frame makes sure you have the strength and vibration-dampening ride to go with your beauty. Style and strength to take on whatever life throws at you.





















Weight capacity— 300 lb

## Specifications

Max user weight	265 lb standard 300 lb heavy duty	Center of gravity adjustment	2.25"
Frame features	Titanium dual-tube rigid, TiFit, limited adjustability	Front seat width (seat taper)	0", 1.5", 2.5"
Total chair width*	23"	Precision lock folding seat back	Adj. height: 9.5" – 22" Adj. angle: 80° – 101°
Seat width	10" – 20"	Fixed seat back	Fixed height: 8" – 20" Fixed angle: 80° – 101°
Seat depth	10" – 20"	Ergonomic seat	4" – 10"
Front seat height	16" – 21"	Transit options	Transport bracket
Rear seat height	12.5" – 20"	Transport weight*	9 lb
Front frame angle	70°, 75°, 80°, 85°, or 90°	Warranty	Five years on frame One year on TiLite components
Camber	0°, 2°, 4°, 6°, 8°, or 12°		








\* Weights and widths based on 16" x 16" chair with 0° camber, minimum rear wheel spacing and all standard features. Rear wheels are not included in transport weight. For reference, a pair of standard 24" rear wheels weighs 7.2 lb.

## Frame colors

High gloss paint			Candy paint		Metallic paint			
								
TiLite red	Neon pink	High gloss black	Ruby red	Azure blue	Electric plum	Tangerine metallic	Acid green metallic	Ocean blue metallic
Matte textured paint			Matte paint		Pearlescent paint		Natural titanium	
								
Flame orange	OD green	Hyper blue	Titanium grey	Black lava	Midnight blue pearl	White pearl	Satin	Polished

## Component colors

# Anodized

		
Red	Orange	Gold
		
Green	Blue	Purple
		
Black		



# TiLite Aero X

The Aero X features an innovative modular design with clean modern lines and configurability. Featuring multiple configurations that include heavy duty, hemi-height, and amputee options. Add superior rolling dynamics and full adjustability, and you get a chair that excels in both performance and price.



















## Specifications

Max user weight	265 lb standard 350 lb heavy duty	Front frame angle	Fixed front: 70°, 80°, or 85° Swing-away: 70°, 80°, or 85°
Frame features	Aluminum box folding frame, modular design, full adjustability	Camber	0°, 2°, 4°, or 6°
Total chair width folded*	10.5"	Center of gravity adjustment	6"
Total chair width*	22.5"	Fixed seat back	Fixed height: 11" – 20"
Seat width	12" – 22"	Transit options	Transit tie-down, transport bracket
Seat depth	12" – 20"	Transport weight*	Fixed front: 15.6 lb Swing-away: 18 lb
Front seat height	13.5" – 21.5"	Warranty	Five years on frame One year on TiLite components
Rear seat height	13" – 20"		

\* Weights and widths based on 16" x 16" chair with 0° camber, minimum rear wheel spacing and all standard features. Rear wheels are not included in transport weight. For reference, a pair of standard 24" rear wheels weighs 7.2 lb.

## Frame colors

High gloss paint			Candy paint		Metallic paint			
								
TiLite red	Neon pink	High gloss black	Ruby red	Azure blue	Electric plum	Tangerine metallic	Acid green metallic	Ocean blue metallic
Matte textured paint			Matte paint		Pearlescent paint			
								
Flame orange	OD green	Hyper blue	Titanium grey	Black lava	Midnight blue pearl	White pearl		





# TiLite 2GX

The TiLite 2GX combines the advantages of a modular folding frame with the high strength and vibration-dampening characteristics of titanium. Add in the Tru-Fit full adjustability that meets your needs and you are ready to roll in a chair designed to empower you.



## Specifications

Max user weight	265 lb standard 350 lb heavy duty	Front frame angle	Fixed front: 70°, 80°, or 85° Swing-away: 70°, 80°, or 85°
Frame features	Titanium box folding frame, modular design, full adjustability	Camber	0°, 2°, 4°, or 6°
Total chair width folded*	10.75"	Center of gravity adjustment	6"
Total chair width*	23"	Fixed seat back	Fixed height: 11" – 20"
Seat width	18" – 22"	Transit options	Transport bracket
Seat depth	16" – 20"	Transport weight*	Fixed front: 13.1 lb Swing-away: 15.5 lb
Front seat height	13.5" – 21"	Warranty	Five years on frame One year on TiLite components
Rear seat height	13" – 20"		

\* Weights and widths based on 16" x 16" chair with 0° camber, minimum rear wheel spacing and all standard features. Rear wheels are not included in transport weight. For reference, a pair of standard 24" rear wheels weighs 7.2 lb.

## Frame colors

High gloss paint			Candy paint			Metallic paint		
TiLite red	Neon pink	High gloss black	Ruby red	Azure blue	Electric plum	Tangerine metallic	Acid green metallic	Ocean blue metallic
Matte textured paint			Matte paint		Pearlescent paint	Natural titanium		
Flame orange	OD green	Hyper blue	Titanium grey	Black lava	Midnight blue pearl	White pearl	Satin	Polished





# SmartDrive MX2+

SmartDrive’s intuitive controls add a power assist for a variety of manual wheelchair scenarios like cruising up hills or over thick carpets. Its lightweight design attaches quickly and its unique OmniWheel design moves with you as you pivot, turn, wheelie, and navigate curbs in style. SmartDrive gives you freedom to do more.



Weight capacity— 331 lb

## Specifications

Max user weight	30 – 331 lb
Product dimensions	Length: 15.3” Width: 5.7” Height: 9.5”
Product weight	13.2 lb drive unit
Driving range*	12.3 mi level/no resistance
Motor power	250 W brushless DC
Maximum speed	5.5 mph
Operating temp	-13°F – 122°F
Chair type compatibility	Rigid, folding, tilt-in-space, one-arm drive, and standing

Drive wheel diameter	22”, 24”, 25”, 26”
Warranty	2-Year limited warranty
Bluetooth® LE	FCC ID   QQBT113 IC ID     5123A-BGTBLE113


*\*Performance measurements will vary based on user and chair characteristics, driving and battery conditions.*

*The SmartDrive has been tested and conforms to all applicable requirements of ANSI/RESNA Standards for Wheelchairs - Volume 1 and Volume 2 and EN 12184.*

*\*\*A wheelchair’s specifications provided by its manufacturer may be slightly affected by the addition of the SmartDrive.\*\**



# Cushions & back supports



### ROHO Contour Select

No weight limit when properly sized to cushion and individual.



ROHO® Hybrid Elite® SR® Single Compartment cushion\*



ROHO Hybrid Elite Dual Compartment cushion\*



ROHO CONTOUR SELECT® cushion

Ideal for individuals who

May require accommodation of pelvic asymmetry, up to approx. 2"  
May require pelvic and lower extremity positioning and a more stable environment for transfers

May require correction of pelvic asymmetry (medial lateral), up to approx. 2"  
May require a more stable environment for transfers

May require accommodation of pelvic asymmetry, up to approx. 1"  
May need correction of pelvic asymmetry (medial lateral AND anterior posterior)  
May require different positioning of the left and right femur

Dimensions

Width: 13.75" – 23.75"  
Depth: 14.75" – 22.75"  
Height: 4.25" without load

Width: 13.75" – 23.75"  
Depth: 14.75" – 22.75"  
Height: 4.25" without load

Width: 14.5" – 24.25"  
Depth: 15.25" – 22"  
Height: 2.25" – 4.25"

Weight limit

Available up to 500 lb

Available up to 500 lb

No weight limit on properly sized cushions

*\*Hybrid Elite images do not show the cushion cover, which is required for use.*



ROHO QUADTRO SELECT® HIGH PROFILE® cushion



ROHO QUADTRO SELECT® MID PROFILE® cushion



ROHO ENHANCER® cushion

Ideal for individuals who

May be at risk for or currently have skin/soft tissue breakdown and/or any category/stage pressure injury  
May require accommodation or correction of medial lateral and/or anterior posterior pelvic asymmetries, up to approx. 3"  
May require different positioning of the left and right femur

May be at risk for or currently have skin/soft tissue breakdown and/or any category/stage pressure injury  
May require accommodation or correction of medial lateral and/or anterior posterior pelvic asymmetries, up to approx. 2"  
May require different positioning of the left and right femur

May require accommodation of pelvic asymmetry, up to approx. 1"  
Require positioning of the pelvis and lower extremities  
May need increased lateral support  
May need to partially offload the ischial tuberosities

Dimensions

Width: 11" – 26.25"  
Depth: 11.75" – 23.75"  
Height: 4.25"

Width: 11" – 26.25"  
Depth: 11.75" – 23.75"  
Height: 3.25"

Width: 14.75" – 20"  
Depth: 14.75" – 20"  
Height: 2.25" – 4.25"




Weight limit

No weight limit on properly sized cushions



No weight limit on properly sized cushions



No weight limit on properly sized cushions



			
	<b>ROHO HIGH PROFILE Dual Compartment cushion</b>	<b>ROHO HIGH PROFILE Single Compartment cushion</b>	<b>ROHO MID PROFILE Single Compartment cushion</b>
<b>Ideal for individuals who</b>	May be at risk for/or currently have skin/soft tissue breakdown and/or any category/stage pressure injury	May be at risk for/or currently have skin/soft tissue breakdown and/or any category/stage pressure injury	May be at risk for/or currently have skin/soft tissue breakdown and/or any category/stage pressure injury
	May require accommodation of pelvic asymmetry, up to approx. 3"	May require accommodation of pelvic asymmetry, up to approx. 3"	May require accommodation of pelvic asymmetry, up to approx. 2"
<b>Dimensions</b>	Width: 14.75" – 20"	Width: 11.5" – 26.75"	Width: 11.5" – 26.75"
	Depth: 14.75" – 21.75"	Depth: 11.5" – 21.75"	Depth: 11.5" – 21.75"
<b>Weight limit</b>	Height: 4.25"	Height: 4.25"	Height: 3.25"
	No weight limit on properly sized cushions	No weight limit on properly sized cushions	No weight limit on properly sized cushions

			
	<b>Comfort Saddle 7 Series cushion</b>	<b>Comfort M2 cushion</b>	<b>Comfort M2 cushion with ATI</b>
<b>Ideal for individuals who</b>	May need postural alignment	May need coccyx relief	May need coccyx relief
	May need relief in the coccyx area	May need postural alignment	May need postural alignment
<b>Dimensions</b>	May require minimal to moderate skin protection	May need support to the greater trochanters	May need support to the greater trochanters
	Requires/wants cushion to strap to chairs		May need dynamic positioning
<b>Weight limit</b>		Width: 18" – 24"	Width: 18" – 24"
		Depth: 14" – 24"	Depth: 14" – 24"
<b>Weight limit</b>		Height: 4.25" for Zero elevation, 3.75" for Anti-Thrust	Height: 4.25" for Zero elevation, 3.75" for Anti-Thrust
		18" wide = 300 lb	18" wide = 300 lb
<b>Weight limit</b>		20" wide = 350 lb	20" wide = 350 lb
		22" wide or > = 400 lb	22" wide or > = 400 lb

			
	<b>Comfort Embrace cushion</b>	<b>Comfort Embrace cushion with ATI</b>	<b>Comfort Ascent cushion</b>
<b>Ideal for individuals who</b>	May need postural alignment	May need postural alignment	May need support for the greater trochanters
	May need support to the greater trochanters	May need support to the greater trochanters	May require mild-moderate skin protection
<b>Dimensions</b>	May need dynamic positioning	May need dynamic positioning	
	Width: 18" – 24"	Width: 18" – 24"	Width: 16" – 24"
<b>Weight limit</b>	Depth: 14" – 24"	Depth: 14" – 24"	Depth: 16" – 20"
	Height: 4" for Zero elevation, 3.75" for Anti-Thrust	Height: 4" for Zero elevation, 3.75" for Anti-Thrust	Height: 4"
<b>Weight limit</b>	18" wide = 300 lb	18" wide = 300 lb	18" wide = 300 lb
	20" wide = 350 lb	20" wide = 350 lb	20" wide = 350 lb
<b>Weight limit</b>	22" wide or > = 400 lb	22" wide or > = 400 lb	22" wide or > = 400 lb

			
	<b>Comfort ACTA-BACK® back support</b>	<b>Comfort ACTA-BACK Deep back support</b>	<b>Comfort ACTA-BACK Contour back support</b>
<b>Ideal for individuals who</b>	May require a basic amount of postural support	May require a significant amount of postural support	
	Those who are more independent with transfers	Those who are more dependent with transfers	May require lower thoracic trunk support
<b>Dimensions</b>	Those who have good trunk control	May need the maximum amount of trunk support	
	Width: 20" – 24"	Width: 20" – 24"	Width: 20" – 24"
<b>Weight limit</b>	Length: 14" – 22"	Length: 14" – 22"	Length: 16" – 22"
	20" wide = 300 lb	20" wide = 300 lb	20" wide = 300 lb
<b>Weight limit</b>	22" wide = 350 lb	22" wide = 350 lb	22" wide = 350 lb
	24" wide = 400 lb	24" wide = 400 lb	24" wide = 400 lb





Comfort ACTA-EMBRACE®  
back support



Comfort ACTA-RELIEF®  
back support

Ideal for individuals who	Complex user that truly needs a custom fit	
	May require a significant amount of postural support	Those who have moderate to severe kyphosis
Dimensions	Width: 20" – 24" Length: 14" – 22"	Width: 20" Length: 10" – 22"
	20" wide = 300 lb 22" wide = 350 lb 24" wide = 400 lb	20" wide = 300 lb
Weight limit		

©2022 Permobil. The following are trademarks and registered trademarks of Permobil: ROHO®, Hybrid Elite®, CONTOUR SELECT®, QUADTRO SELECT®, HIGH PROFILE®, MID PROFILE®, ENHANCER®, ACTA-BACK®, ACTA-EMBRACE®, ACTA-RELIEF®.

Permobil has a policy of continual product improvement and reserves the right to amend this document. The current version of this document is available at permobil.com.



Not seeing your size?

Contact your Permobil rep for more details.

## Other offerings



Bodypoint belts



Rifton E-Pacer



Rifton TRAM

Product overview	Bodypoint padded and non-padded hip belts come in extra-large 81" lengths and are comfortable, simple-to-use, and easy to clean.	The E-Pacer has a sit-to-stand lift functionality that enables gait training with large or highly dependent clients. Its strong and stable frame can accommodate users up to 6' 5".	Capable of lifting and transferring clients of varying abilities and body types, the TRAM is easy to maneuver in tight spaces and simple to operate. Can also function as a gait training device.
	Weight limit	—	350 lb




## Tröst Plus bed

With a weight capacity of up to 600 lb, the Tröst Plus adjustable bed is ideal for heavy weight clients. It offers a low height of 6.75" and a high height of 30" to meet all your care-giving and client needs. Tröst has a 12-button hand pendant and customized comfort is an easy adjustment. Available in 39" or 42" widths and extendable to 84" in length.



## Hälsa Plus bed

Hälsa Plus is an electric, adjustable home care bed with a weight capacity of 500 lb. It offers countless head and foot positioning options with the 10-button hand pendant and is easy to disassemble for transportation and storage. Available in four width options: 36", 39", 42" or 48".



## Resten Select bed

The new Resten Select adjustable bedframe can accommodate almost any position to help with better sleep, improved circulation, and decreased back pain. Key features include programmable massage, memory settings, USB charging ports built-in, and a wireless controller. Queen size weight capacity is up to 600 lb.





## Serene Air

Serene Air is a dual mode therapeutic surface with Alternating Pressure Therapy and Pulsation options that assist in wound healing and prevention. It offers 120 lpm of Low Air Loss for effective microclimate control and supports clients weighing up to 920 lb. Additional features include seat inflation for the Fowler position, a quilted top cover, and a choice of APT operating cycle times. Available in three widths: 36", 42" or 48".




## Domus Edge

Domus Edge Alternating Pressure Therapy system takes powered therapeutic surfaces to the next level with the added safety feature of a 3" wide transfer edge. This helps enable safe and easy transfers while also reducing falls and entrapment for clients weighing up to 450 lb. System has seat inflation for the Fowler position, a quilted top cover, and lower static cells that prevent clients from bottoming out. Available in 36" and 42" widths.



## Adapt 2

Adapt 2 Series modular mattress systems are an innovative way to reduce the effects of pressure, friction, and shear. Up to four foam inserts or ROHO mattress sections can be added to the surface based on the individual client's needs. Available in three widths: 36", 39" or 42" and has a weight capacity of 350 lb.



## Align Bariatric

This unique multi-zone surface with Shear Reduction Technology and firm borders is an ideal sleep surface for clients weighing 600-800 lb. There are differing immersion levels that match the contours of the body, and the washable Recovery 5 cover provides incontinent protection. Available in six widths from 36" – 60".