

Mattresses & therapeutic surfaces matrix

Features																
	Non-powered pressure redistrib.	Shear Reduction Technology (SRT)	Insert ROHO mattress sections	Gel swirl foam	Firm side borders	Raised bolster custom option	Special foot section	Alternating Pressure Therapy (APT)	Low Air Loss	Pulsation therapy	Continuous Lateral Rotation Technology	Recovery 5 or Recovery 5 Ultra cover	Fluid resistant cover	Plus/ bariatric sizes	Weight limit (lb)	Warranty
Mattresses																
Adapt	✓	✓	✓	✓	✓	✓	✓					✓		✓	350-500	24 months
Allay	✓	✓		✓	✓	✓						✓		✓	350-500	24 months
Align & Align Bari.	✓	✓			✓	✓	✓					✓		✓	350-800	24 months
Sömn	✓				✓	✓							✓ Startex	✓	350-500	24 months
Koppla	✓				✓								✓ Startex	✓	350-500	24 months
Optima 2	✓													✓	350-500	24 months
Symmetry	✓				✓		✓						✓ stretch	✓	350-500	24 months
Resten Luxe				✓										✓	350-600	24 months
Therapeutic Surfaces - Powered (Active)																
Domus 4 w/ foam edge					✓			✓	✓ 8 L/min				✓ quilted	✓	450	36 months
Domus 4								✓	✓ 8 L/min				✓ quilted		450	36 months
Domus Auto								✓	✓ 8 L/min				✓ quilted		450	36 months
Serene Air w/ foam edge					✓			✓	✓ 120 L/min	✓			✓ quilted	✓	400-650	36 months
Serene Air								✓	✓ 120 L/min	✓			✓ quilted	✓	400-920	36 months
Pro-Care Optima							✓	✓	✓ 13 L/min				✓ stretch		550	36 months
Pro-Care Turn							✓	✓	✓ 13 L/min		✓		✓ stretch		395	36 months
Therapeutic Surfaces - Non-Powered (Reactive)																
Atlas Hybrid	✓			✓	✓		✓						✓ stretch	✓	500-700	36 months
Atlas Auto Comfort	✓			✓	✓	✓							✓ stretch	✓	500-700	36 months
Vicair 415	✓												✓		550	lifetime on SmartCells
Therapeutic Surfaces - Modular Non-Powered																
Adapt 2 Series	✓	✓	✓	✓	✓		✓					✓		✓	350	24 months
Acute	✓		✓	✓			✓					✓			350	24 months
Mattress Overlays																
ROHO Prodigy	✓														300	12 months
ROHO SOFFLEX 2	✓														330	18 months
ROHO Dry Floatation	✓												✓	✓	unlimited	24 months
Leveling pads	✓		✓									✓			250-1,000	24 months