


SCOUT 4



1. Speed Dial- BC31001

Speed Dial Knob- LRM412103-11

2. Throttle Assembly Old [2]

Throttle Assembly New- S12-069-00300

3. Battery Gauge (4 LEDs)- SC31004

Battery Gauge (8 LEDs)- C11-076-00100 [1]

4. Horn Button- SC31034

5. Key Switch Assembly- BC31006

Key- BCKEY

6. Upper Tiller Cover- SC31008

7. Lower Tiller Cover- SC31009

8. Basket Bracket- BC35051

9. Basket- S35003-P

10. Handgrip Foam- SC31002

Handgrip Plastic- SC31054

11. Tiller- SC31001

12. Tiller Swivel Assembly- SC31010

13. Tiller Lock- C08-069-01000

14. Tiller Bearing Set- BC31021

15. "L" Steering Stem- SC31016

16. Steering Fork Assembly- SC31015

17. Steering Axle (left or right)- SC31021

18. Connecting Rod- BC350103

19. Bumper Bracket- SC31022

20. Front Bumper- SC31028

21. Front Wheel Black- SC31017

Deluxe Front Wheel- SC31017-D

22. Front Shroud- SC31013

23. Footrest Shroud- SC31024

24. Battery Box Shroud- SC31027

25. Red Decoration- SC31029

Deluxe Red Decoration- SC31029-D

Blue Decoration- SC31030

Deluxe Blue Decoration- SC31030-D

26. Front Frame- SC31014

27. Black Seat- SC35045

Deluxe Two-Tone Seat- SC31035

Deluxe Back Cover- SC31035-BC

Deluxe Seat Cover- SC31035-SC
28. Left Armrest- BC35040L

Right Armrest- BC35040R

Arm Pad- BC953014

29. Armrest Knob- S31052

30. Seat Plate- SC31023

31. Seat Post- SC31018 [4]

32. Seat Post Bolt- C07-071-00101

33. Battery Box & Batteries (12V 12ah)- SC31007

Battery Box Only (12V 12ah)- SC31006

Battery 12V 12ah (2 needed)- LRM402109

EXT Model

Battery Box & Batteries (12V 20ah)- SC31050

Battery Box Only (12V 20ah)- SC31051

Battery 12V 20ah (2 needed)- LRM412123

34. Charger Socket & Harness- SC31031

35. 20 Amp Circuit Breaker- BC20ACB

36. Battery Contact- C12-060-00300

37. Disconnect Handle Assembly- SC31019

38. Motor/Brake/Transaxle Assembly- SC31053

Transaxle- BC350112

Motor- SC31036

Electro Magnetic Brake- SC31037

39. Motor Wire Harness- SC31020

40. Rear Wheel Black- SC31003

Deluxe Rear Wheel- SC31003-D

41. Controller Curtis- SC31052-4

Controller Dynamic- C10-069-00100-4 [3]

42. Wire Harness Complete-

Curtis- SC31026

Dynamic- DST-018 [3]

43. Contact Block Set w/Motor Harness- SC31033

44. Rear Shroud- SC31011

45. Rear Frame- SC31012

46. Anti-Tipper Wheel- S31026

47. Charger- LRM302111

48. Front Wheel Bearing- S350100-B

49. Hub Cap- SC31025

50. Axle Bolt Assembly- SC31049

51. Battery Box Hold Down- SC31047

[1] After Serial # VP035868
[2] Stock is gone use new style.
[3] After Serial # VP142824
[4] No longer available in chrome once stock is gone.